

PHS-25CW/3BW Benchtop pH Meter



Features:

- Economical benchtop pH meter is equipped with a white backlit LCD display.
- 2 points push-button calibration with auto-buffer recognition.
- Automatic electrode slope display helps user decide whether to replace sensor.
- Automatic Temperature Compensation ensures accurate measuring result.

Meter Includes:

- pH electrode
- Temperature probe (Only for PHS-3BW meter)
- pH buffer pouches (pH4.01/7.00/10.01)
- DC9V power adapter

Specifications:

Model	PHS-25CW	PHS-3BW
pH Range	0.00~14.00pH	0.00~14.00pH
pH Accuracy	±0.05pH	±0.01pH
Resolution	0.01pH	0.01pH
Calibration Points	2 points	2 points
Calibration Solutions	USA (pH4.01/7.00/10.01)	USA (pH4.01/7.00/10.01)
Temperature Compensation	0~100°C, Manual	0~100°C, Manual or Automatic
mV Range	-1999~1999mV	-1999~1999mV
mV Accuracy	±10mV	±1mV
Resolution	1mV	1mV
Temperature Range	—	0~100°C
Temperature Accuracy	—	±1°C
Resolution	—	1°C
Hold Function	—	—
Connector	BNC	BNC
Display	LCD	LCD
Power Requirements	DC9V/400mA power adapter	DC9V/400mA power adapter
Dimensions	210(L)×205(W)×65(H)mm	210(L)×205(W)×65(H)mm
Weight	1.5Kg	1.5Kg