

## Thermo Scientific Orion Star A211 pH Benchtop Meter

The Thermo Scientific Orion Star A211 pH Benchtop Meter works for any lab or application. Its durability and dependability make it ideal for every sample, every day.



The Thermo Scientific Orion Star A211 pH Benchtop Meter is the ideal choice for any lab. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV and temperature along with electrode status, time, date, sample ID, user ID and calibration points. For advanced users, features such as stability and averaging choices provide additional options.

### Features and Benefits

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading – AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix calibration errors without a complete recalibration with calibration editing
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Mix samples without a stir plate with direct control of the Orion Star stirrer probe (096019, sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty

## Thermo Scientific Orion Star A211 pH Benchtop Meter

### Specifications



<b>pH</b>	<b>Range</b>	-2.000 to 20.000
	<b>Resolution</b>	0.1, 0.01, 0.001
	<b>Relative Accuracy</b>	±0.002
	<b>Calibration Points</b>	Up to 5
	<b>Calibration Editing</b>	Yes
<b>mV/RmV</b>	<b>Range</b>	±2000.0 mV
	<b>Resolution</b>	0.1
	<b>Relative Accuracy</b>	±0.2 mV or ±0.05 % of reading whichever is greater
	<b>E<sub>H</sub> ORP Mode</b>	Yes
<b>Temperature</b>	<b>Range</b>	-5 to 105 °C, 23 to 221 °F
	<b>Resolution</b>	0.1
	<b>Relative Accuracy</b>	±0.1
	<b>Offset Calibration</b>	1 point
<b>Datalogging</b>	<b>Memory</b>	2000 with time and date stamp
	<b>Log Function</b>	Manual, ready (includes AUTO-READ™), timed
	<b>Log Edit</b>	Delete individual points, group of points or all
<b>Inputs</b>	<b>pH Electrode</b>	BNC, reference pin
	<b>ATC Probe</b>	8-pin mini-DIN
<b>Output</b>		RS232, USB
<b>Power</b>	<b>AC Adapter</b>	Included – universal, 100-240 VAC
	<b>Battery Power</b>	Optional – 4 AAs
	<b>Battery Life</b>	800 hrs

### Ordering Information

Cat. No.	Description
<b>STARA2110</b>	<b>Orion STAR A211 pH Benchtop Meter</b> - Electrode arm with redesigned holder - Universal power adapter
<b>STARA2115</b>	<b>Orion STAR A211 pH Benchtop Meter Kit</b> - 8302BNUMD Orion ROSS Triode 3-in-1 pH/ATC probe - 810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle) - Electrode arm with redesigned holder - Universal power adapter
<b>Accessories</b>	
<b>096019</b>	Orion Star stirrer probe
<b>STARA-BEA</b>	Benchtop electrode arm for Orion Star A-series meters
<b>STARA-HB</b>	Freestanding base for use with Orion Star A-series benchtop electrode arm
<b>8302BNUMD</b>	Orion ROSS Triode 3-in-1 pH/ATC probe, refillable, glass body
<b>8102BNUWP</b>	Orion ROSS Ultra refillable pH electrode, glass body
<b>927007MD</b>	Orion ATC probe, stainless steel
<b>810199</b>	Orion ROSS All-in-One pH buffer kit
<b>910104</b>	Orion pH 4.01 buffer, 475 mL
<b>910107</b>	Orion pH 7.00 buffer, 475 mL
<b>910110</b>	Orion pH 4.01 buffer, 475 mL

For more information, contact your local Thermo Scientific products dealer or call our customer and technical service experts at 1-800-225-1480 (for the US and Canada) or visit [www.thermoscientific.com/water](http://www.thermoscientific.com/water).

©2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. ROSS and the COIL tradenames are trademarks of Thermo Fisher Scientific Inc. US Patent 6,793,787.



S-STARA211-E-1011-RevB

**Water Analysis Instruments**  
[www.thermoscientific.com/water](http://www.thermoscientific.com/water)

**North America**  
166 Cummings Center  
Beverly, MA 01915 USA  
Toll Free: 1-800-225-1480  
Tel: 1-978-232-6000  
info.water@thermo.com

**Netherlands**  
Tel: (31) 033-2463887  
info.water.uk@thermo.com

**China**  
Tel: (86) 21-68654588  
wai.asia@thermofisher.com

**India**  
Tel: (91) 22-4175-8800  
wai.asia@thermofisher.com

**Singapore**  
Tel: (65) 6778-6876  
wai.asia@thermofisher.com

**Japan**  
Tel: (81) 045-453-9175  
wai.asia@thermofisher.com

**Australia**  
Tel: (613) 9757-4300  
infowaterau@thermofisher.com

**Thermo**  
SCIENTIFIC